

Attachment B - Pricing Sheet
Solicitation # 3000023542
VFMC Wound Vac Rental and Supplies

Line #	Item Description	Estimated Quantity	Unit of Measure	Unit Price	Extended Amount
1	<p>Rental of a Portable Negative Pressure Wound Therapy (NPWT) device for advanced wound healing (*NOTE: This device must be compatible with all of the items being bid on lines #3-#22.) Device being bid must meet all of the following specifications:</p> <ul style="list-style-type: none"> - hold a minimum 300 mL volume disposable canister (to minimize canister changes) - have an on-screen user guide - be portable, to allow for patient mobility and to include a carrying case with a pocket to store excess tubing for on-the-go therapy - include a battery with a minimum 12-hour life (long enough to enable patient to be mobile for a full day) - have an adjustable rate of dressing draw down intensity for increased patient comfort - have a filter system (to minimize wound odor) - include technology (or installed software) with capabilities to continuously monitor, measure and maintain the prescribed negative pressure at the wound site for optimal healing outcomes by adjusting pump output, compensating for wound distance, wound position, exudate characteristics, and patient movement that includes safety alarms for patient safety to alert caregivers if target pressure is not met or therapy is interrupted. This technology must be included in the tubing and pad to efficiently draw exudate away from the wound, independently monitor target pressure at the wound, and to distribute negative pressure to help reduce blockages and false alarms, and to enhance patient comfort <p>Dimensions: Approx. 6.5"W X 3.125"H X 7.5"D Approx weight: 3.20 lbs.</p> <p>Therapy Adjustment Options: Negative pressure: min. 50-200 mmHG Mode of operation: Continuous or Intermittent Intensity: Variable</p> <p>Electrical: Charger/External Supply Input: 100-240 VAC, 50/60 Hz 1 Amp Charger/External Supply Output: 18 VDC, 5 Amp Patient and Enclosure Leakage Current (max.): <100 microamps</p> <p>Storage Conditions: Temperature Range: -4 degrees to 140 degree</p> <p>Operating Conditions: Temperature Range: 50 degrees F to 85 degrees F</p> <p>Display Screen Features must include, but are not limited to: Power, Exit, Select, Enter and Arrow keys Therapy, Utilities and Language options Alarm silence, AC power and battery indicators LED lights</p> <p>***Must have a lock-out feature***</p> <p>Alarm type features must include, but not limited to: Canister full, tubing blocked, canister missing, tubing leaks, dressing leaks, therapy not activated, audible low battery, internal device, and internal memory error.</p> <p>IEC Classification: Medical Equipment Continuous Operation Type B Applied Part Class II Internally Powered Equipment IPX0</p> <p>List contact phone numbers for ordering and scheduling pick-up of devices and when calling to let vendor know when device is being put into use:</p> <p>Ordering: _____</p> <p>Scheduling Pick-up: _____</p> <p>KCI V.A.C. Freedom Therapy Unit, Item #320000.P or equal</p> <p>State canister size (mL): _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Day		

Attachment B - Pricing Sheet
Solicitation # 3000023542
VFMC Wound Vac Rental and Supplies

Line #	Item Description	Estimated Quantity	Unit of Measure	Unit Price	Extended Amount
2	<p>Rental of a Negative Pressure Wound Therapy (NPWT) device that provides three therapies all by one device. (*NOTE: This device must be compatible with all of the items being bid on lines # 3-20 and # 23-31)</p> <p>Device must be one unit that provides all of the following three therapies:</p> <ol style="list-style-type: none"> 1) NPWT, 2) Instillation Therapy (automated delivery of topical wound solutions), and 3) Open Abdomen Negative Pressure Therapy. <p>Device must hold at least two different size disposable canisters: min. 500 mL and min. 1000 mL.</p> <p>Device must have an onscreen user guide/touch screen menu, simple therapeutic programming, and audible and visual alarms.</p> <p>Device must include installed technology (or software) with capabilities to continuously monitor, measure and maintain prescribed negative pressure at the wound site for optimal healing outcomes by adjusting pump output, compensating for wound distance, wound position, exudate characteristics, and patient movement that includes safety alarms to alert caregivers if target pressure is not met or therapy is interrupted. This technology must be included in the tubing and pad to efficiently draw exudate away from the wound, independently monitor target pressure at the wound, and to distribute negative pressure to help reduce blockages and false alarms, and to enhance patient comfort.</p> <p>Device must have all of the following:</p> <ul style="list-style-type: none"> -technology for leak detections (to offer instant feedback to help identify leaks). -technology providing continuous feedback for enhanced negative pressure accuracy, monitoring pressure at the wound site, and adjusts accordingly for optimal healing outcomes. -features to help clinicians identify and troubleshoot negative pressure leaks, provide real-time feedback with audible and visible cues for leak location, and guide clinicians to achieve a proper seal. <p>Device must be capable of containing history reports to allow clinicians to monitor alarm history, patient history, and therapy history.</p> <p>Device must have customized pressure algorithms: continuous negative pressure (from -25 to -200mmHg), technology to maintain minimum negative pressure levels (-25mmHg) between cycles, helping to prevent leaks and fluid accumulation.</p> <p>Dimensions: Approx. 8.55"W X 10.25"H X 7.5"D Approx. weight: 7.4 lbs</p> <p>Electrical: Voltage: 100-240 VAC Frequency: 50Hz-60Hz Power: 60W Battery Life: approx. 6 hours</p> <p>Environmental Conditions: Transport & Storage Temperature Range: -4 degrees F to 140 degrees F Operational Temperature Range: 50 degrees F to 86 degrees F</p> <p>Device must have Touch Screen Features such as, but not limited to: Status Bar, Tabs for: Therapy, History and Utilities, Therapy Unit Configuration and Phase Status, Therapy Settings, Leak Detector, Therapy Phase, Dressing Soak, Therapy Pause/Resume and Start/Stop, Test Cycle, Help, Night Mode, Screen Guard/Settings Lock, Cancel.</p> <p>Must have Other Features such as, but not limited to: Battery charging indicating light, Speaker, SD card slot, USB port stylus. Solution container hanger arm, unit hanger arm, solution container hanger arm lock. Carry handle, Cassette tubing retention groove, Cassette pivot slot. Canister bellows (seals).</p> <p>Device must have Therapy Settings such as, but not limited to: Default settings.</p> <p>Able to set fill assist (off/on: to help user determine the correct instill volume to use during the instill phase of a therapy cycle), start phase, instill volume, soak time (minutes), therapy time, target pressure, intensity (low, medium, high).</p> <p>Low Priority/Alert Conditions: Must display on the touch screen when the therapy unit detects a condition that requires attention and are accompanied by repeating audible tones options such as, but not limited to: Audible pause, battery low, therapy blockage, solution bag/bottle empty, fill assist inactive, internal temperature.</p>	1	Day		

Line #	Item Description	Estimated Quantity	Unit of Measure	Unit Price	Extended Amount
	<p>Medium Priority Alarm Conditions: Must display on the touch screen when the therapy unit detects a condition that requires prompt attention and are accompanied by a repeating audible tone and a flashing screen title options such as, but not limited to: Battery critical, canister full, canister not engaged, therapy leak, pressure deviation, therapy blockage, therapy low pressure, system error, therapy inactive.</p> <p>IEC Classification: Continuous operation Class II or internally powered equipment Type BF applied part IPX1</p> <p>State contact numbers for ordering and scheduling pick-up of devices: Ordering: _____ Scheduling Pick-up: _____</p> <p>KCI V.A.C. Ulta Therapy Unit, Item #ULTDEV01/US.P or equal</p> <p>State the two different canister sizes (mL): Smaller canister (mL): _____ Larger canister (mL): _____</p> <p>Specify brand (and number if applicable): _____</p>				
Line Items #3-20 are Disposables or Accessories to purchase to go with the devices being bid on both lines #1 and #2. Therefore, all items must be compatible with lines #1 and #2.					
3	<p>Disposable, Standard, Sterile Foam Dressing Kit, Size Small (to go with both Line Items #1 and #2): Individually wrapped sterile packaging, for one-time use Included in each kit: *One black foam dressing (dimensions vary according to kit size; can be cut to size) *One, two or five sheet(s) of standard clear, occlusive drape (approx. 30.5x26 cm: number of sheets included depends on kit size; can be cut to size) *One pad with tubing (including clamp and connector which allows two dressings to be connected to one therapy unit) *One disposable ruler Foam is 3" - 5" thick and can be cut to fit. Dressing size: approx. 10 X 7.5 X 3.2 cm Contains at least one sheet of drape Packaged: 10 kits/case</p> <p>KCI V.A.C. GranuFoam Small Dressing, Item #M8275051/10.S or equal</p> <p>State number of kits/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
4	<p>Disposable, Standard, Sterile Foam Dressing Kit, Size Small (to go with both Line Items #1 and #2): Individually wrapped sterile packaging, for one-time use Included in each kit: *One black foam dressing (dimensions vary according to kit size; can be cut to size) *One, two or five sheet(s) of standard clear, occlusive drape (approx. 30.5x26 cm: number of sheets included depends on kit size; can be cut to size) *One pad with tubing (including clamp and connector which allows two dressings to be connected to one therapy unit) *One disposable ruler Foam is 3" -5" thick and can be cut to fit. Dressing size: 10 X 7.5 X 3.2 cm Contains at least one sheet of drape Packaged: 5 kits/case</p> <p>KCI V.A.C. GranuFoam Small Dressing, Item #M8275051/5.S or equal</p> <p>State number of kits/case: _____</p>	1	Case		

Attachment B - Pricing Sheet
Solicitation # 300023542
VFMC Wound Vac Rental and Supplies

Line #	Item Description	Estimated Quantity	Unit of Measure	Unit Price	Extended Amount
	Specify brand (and number of applicable): _____				
5	<p>Disposable, Standard, Sterile Foam Dressing Kit, Size Medium (to go with both lines #1 and #2): Individually wrapped sterile packaging, for one-time use Included in each kit: *One black foam dressing (dimensions vary according to kit size; can be cut to size) *One, two or five sheet(s) of standard clear, occlusive drape (approx. 30.5x26 cm: number of sheets included depends on kit size; can be cut to size) *One pad with tubing (including clamp and connector which allows two dressings to be connected to one therapy unit) *One disposable ruler Foam is 3" - 5" thick and can be cut to fit. Dressing size: approx. 18 X 12.5 X 3.2 cm Contains at least two sheets of drape Packaged: 10 kits/case</p> <p>KCI V.A.C. GranuFoam Medium Dressing, Item #M8275052/10.S or equal</p> <p>State number of kits/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
6	<p>Disposable, Standard, Sterile Foam Dressing Kit, Size Medium(to go with both lines #1 and #2): Individually wrapped sterile packaging, for one-time use Included in each kit: *One black foam dressing (dimensions vary according to kit size; can be cut to size) *One, two or five sheet(s) of standard clear, occlusive drape (approx. 30.5x26 cm: number of sheets included depends on kit size; can be cut to size) *One pad with tubing (including clamp and connector which allows two dressings to be connected to one therapy unit) *One disposable ruler Foam is 3" - 5" thick and can be cut to fit. Dressing size: approx. 18 X 12.5 X 3.2 cm Contains at least two sheets of drape Packaged: 5 kits/case</p> <p>KCI V.A.C. GranuFoam Medium Dressing, Item #M8275052/5.S or equal</p> <p>State number of kits/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
7	<p>Disposable, Standard, Sterile Foam Dressing Kit, Size Large (to go with both line items #1 and #2): Individually wrapped sterile packaging, for one-time use Included in each kit: *One black foam dressing (dimensions vary according to kit size; can be cut to size) *One, two or five sheet(s) of standard clear, occlusive drape (approx. 30.5x26 cm: number of sheets included depends on kit size; can be cut to size) *One pad with tubing (including clamp and connector which allows two dressings to be connected to one therapy unit) *One disposable ruler Foam is 3" - 5" thick and can be cut to fit. Dressing size: approx. 26 X 15 X 3.2 cm Contains at least two sheets of drape Packaged: 10 kits/case</p> <p>KCI V.A.C. GranuFoam Large Dressing, Item #M8275053/10.S or equal or equal</p> <p>State number of kits/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
8	<p>Disposable, Standard, Sterile Foam Dressing Kit, Size Large (to go with lines #1 and #2): Individually wrapped sterile packaging, for one-time use Included in each kit: *One black foam dressing (dimensions vary according to kit size; can be cut to size)</p>				

Attachment B - Pricing Sheet
Solicitation # 3000023542
VFMC Wound Vac Rental and Supplies

Line #	Item Description	Estimated Quantity	Unit of Measure	Unit Price	Extended Amount
	<p>*One, two or five sheet(s) of standard clear, occlusive drape (approx. 30.5x26 cm: number of sheets included depends on kit size; can be cut to size) *One pad with tubing (including clamp and connector which allows two dressings to be connected to one therapy unit) *One disposable ruler Foam is 3" - 5" thick and can be cut to fit. Dressing size: approx. 26 X 15 X 3.2 cm Contains at least two sheets of drape Packaged: 5 kits/case</p> <p>KCI V.A.C. GranuFoam Large Dressing, Item #M8275053/5.S or equal</p> <p>State number of kits/ case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
9	<p>Disposable, Standard, Sterile Foam Dressing Kit, Size X-Large (to go with both lines #1 and #2): Individually wrapped sterile packaging, for one-time use Included in each kit: *One black foam dressing (dimensions vary according to kit size; can be cut to size) *One, two or five sheet(s) of standard clear, occlusive drape (approx. 30.5x26 cm: number of sheets included depends on kit size; can be cut to size) *One pad with tubing (including clamp and connector which allows two dressings to be connected to one therapy unit) *One disposable ruler Foam is 3" - 5" thick and can be cut to fit. Dressing size: approx. 60 X 30 X 1.5 cm Contains at least 5 sheets of drape Packaged: 5 kits/case</p> <p>KCI V.A.C. GranuFoam Extra Large Dressing, Item #M8275065/5.S or equal</p> <p>State number of kits/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
10	<p>Disposable, Sterile Antimicrobial Foam Dressing Kit, Size Small(to go with both lines #1 and #2): Individually wrapped sterile packaging, for one time use,</p> <p>Foam dressing is impregnated with ionic silver (used to treat already infected wounds or wounds susceptible to infection). Included in each kit: *One silver black foam dressing (dimensions vary according to kit size; can be cut to size) *One or two sheet(s) of standard clear, occlusive drape (approx. 30.5X26 cm: number of sheets included depends on kit size; can be cut to size) *One pad with tubing (including clamp and connector which allows two dressings to be connected to one therapy unit) *One disposable ruler Foam is 3" - 5" thick and can be cut to fit Dressing size: approx. 10 X 7.5 X 3.2 cm Contains at least 1 sheet of drape Packaged: 10 kits/case</p> <p>KCI V.A.C. GranuFoam Silver Small Dressing, Item #M8275098/10.S or equal</p> <p>State number of kits/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
11	<p>Disposable, Sterile Antimicrobial Foam Dressing Kit, Size Small(to go with both lines #1 and #2): Individually wrapped sterile packaging, for one time use,</p> <p>Foam dressing is impregnated with ionic silver (used to treat already infected wounds or wounds susceptible to infection). Included in each kit:</p>				

Attachment B - Pricing Sheet
Solicitation # 3000023542
VFMC Wound Vac Rental and Supplies

Line #	Item Description	Estimated Quantity	Unit of Measure	Unit Price	Extended Amount
	<p>*One silver black foam dressing (dimensions vary according to kit size; can be cut to size) *One or two sheet(s) of standard clear, occlusive drape (approx. 30.5X26 cm: number of sheets included depends on kit size; can be cut to size) *One pad with tubing (including clamp and connector which allows two dressings to be connected to one therapy unit)</p> <p>Foam is 3" - 5" thick and can be cut to fit Dressing size: approx. 10 X 7.5 X 3.2 cm Contains at least one sheet of drape Packaged: 5 kits/case</p> <p>KCI V.A.C. GranuFoam Silver Small Dressing, Item #M8275098/5.5 or equal</p> <p>State number of kits/case: _____ Specify brand (and number if applicable): _____</p>	1	Case		
12	<p>Disposable, Sterile Antimicrobial Foam Dressing Kit, Size Medium (to go with both lines #1 and #2): Individually wrapped sterile packaging, for one time use,</p> <p>Foam dressing is impregnated with ionic silver (used to treat already infected wounds or wounds susceptible to infection). Included in each kit: *One silver black foam dressing (dimensions vary according to kit size; can be cut to size) *One or two sheet(s) of standard clear, occlusive drape (approx. 30.5X26 cm: number of sheets included depends on kit size; can be cut to size) *One pad with tubing (including clamp and connector which allows two dressings to be connected to one therapy unit) *One disposable ruler</p> <p>Foam is 3" - 5" thick and can be cut to fit Dressing size: approx. 18 X 12.5 X 3.2 cm Contains at least two sheets of drape Packaged: 10 kits/case</p> <p>KCI V.A.C. GranuFoam Silver Medium Dressing, Item #M8275096/10.5 or equal</p> <p>State number of kits/case: _____ Specify brand (and number if applicable): _____</p>	1	Case		
13	<p>Disposable, Sterile Antimicrobial Foam Dressing Kit, Size Medium (to go with both lines #1 and #2): Individually wrapped sterile packaging, for one time use,</p> <p>Foam dressing is impregnated with ionic silver (used to treat already infected wounds or wounds susceptible to infection). Included in each kit: *One silver black foam dressing (dimensions vary according to kit size; can be cut to size) *One or two sheet(s) of standard clear, occlusive drape (approx. 30.5X26 cm: number of sheets included depends on kit size; can be cut to size) *One pad with tubing (including clamp and connector which allows two dressings to be connected to one therapy unit) *One disposable ruler</p> <p>Foam is 3" - 5" thick and can be cut to fit Dressing size: approx. 18 X 12.5 X 3.2 cm Contains at least two sheets of drape Packaged: 5 kits/case</p> <p>KCI V.A.C. GranuFoam Silver Medium Dressing, Item #M8275096/5.5 or equal</p> <p>State number of kits/case: _____ Specify brand (and number if applicable): _____</p>	1	Case		
14	<p>Disposable, Sterile Antimicrobial Foam Dressing Kit, Size Large (to go with both lines #1 and #2): Individually wrapped sterile packaging, for one time use,</p> <p>Foam dressing is impregnated with ionic silver (used to treat already infected wounds or wounds susceptible to infection). Included in each kit: *One silver black foam dressing (dimensions vary according to kit size; can be cut to size)</p>				

Attachment B: Pricing Sheet
for Solicitation # 3000023542
VFMC Wound Vac Rental and Supplies

Line #	Item Description	Estimated Quantity	Unit of Measure	Unit Price	Extended Amount
	<p>*One or two sheet(s) of standard clear, occlusive drape (approx. 30.5X26 cm: number of sheets included depends on kit size; can be cut to size) *One pad with tubing (including clamp and connector which allows two dressings to be connected to one therapy unit) *One disposable ruler Foam is 3" - 5" thick and can be cut to fit Dressing size: approx. 26 X 15 X 3.2 cm Contains at least two sheets of drape Packaged: 10 kits/case</p> <p>KCI V.A.C. GranuFoam Silver Large Dressing, Item #M8275099/10.S or equal</p> <p>State number of kits/case: _____ Specify brand (and number if applicable): _____</p>	1	Case		
15	<p>Disposable, Sterile Antimicrobial Foam Dressing Kit, Size Large (to go with both lines #1 and #2): Individually wrapped sterile packaging, for one time use,</p> <p>Foam dressing is impregnated with ionic silver (used to treat already infected wounds or wounds susceptible to infection). Included in each kit: *One silver black foam dressing (dimensions vary according to kit size; can be cut to size) *One or two sheet(s) of standard clear, occlusive drape (approx. 30.5X26 cm: number of sheets included depends on kit size; can be cut to size) *One pad with tubing (including clamp and connector which allows two dressings to be connected to one therapy unit) *One disposable ruler Foam is 3" - 5" thick and can be cut to fit Dressing size: approx. 26 X 15 X 3.2 cm Contains at least two sheets of drape Packaged: 5 kits/case</p> <p>KCI V.A.C. GranuFoam Silver Large Dressing, Item #M8275099/5.S or equal</p> <p>State number of kits/case: _____ Specify brand (and number if applicable): _____</p>	1	Case		
16	<p>Disposable Y connector (to go with both lines #1 and #2): Individually wrapped sterile packaging, for one-time use. Allows two or more dressings to be connected to one therapy unit. Packaged: 10/case</p> <p>KCI T.R.A.C. Y-Connector, Item #M6275066/10.S or equal</p> <p>State number of Y connectors/case: _____ Specify brand (and number if applicable): _____</p>	1	Case		
17	<p>Disposable Y connector (to go with both lines #1 and #2): Individually wrapped sterile packaging, for one-time use. Allows two or more dressings to be connected to one therapy unit. Packaged: 5/case</p> <p>KCI T.R.A.C. Y-Connector, Item #M6275066/5.S or equal</p> <p>State number of Y connectors/case: _____ Specify brand (and number if applicable): _____</p>	1	Case		
18	<p>Disposable, Sterile Tubing Cap (to go with both lines #1 and #2): Individually wrapped sterile packaging, for one-time use. Used to create a seal at the end of canister tubing when not in use . Packaged: 10/case</p> <p>KCI T.R.A.C. Tubing Cap, Item #M6275069/10.S or equal</p> <p>State number of caps/case: _____ Specify brand (and number if applicable): _____</p>	1	Case		

Attachment B: Pricing Sheet
for Solicitation # 3000023542
VFMC Wound Vac Rental and Supplies

Line #	Item Description	Estimated Quantity	Unit of Measure	Unit Price	Extended Amount
19	<p>Disposable, Sterile Tubing Cap (to go with both lines #1 and #2): Individually wrapped sterile packaging, for one-time use. Used to create a seal at the end of canister tubing when not in use . Packaged: 5/case</p> <p>KCI T.R.A.C. Tubing Cap, Item #M6275059/5.S or equal</p> <p>State number of caps/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
20	<p>Disposable, Sterile Drape (to go with both line items #1 and #2): Individually wrapped sterile packaging wrapped sterile packaging, for one-time use. Occlusive clear sheet. Covers the foam dressing(s). Can be cut to size. Standard size: approx. 30.5 X 26 cm Packaged: 10/case</p> <p>KCI V.A.C. Drape, Item #M6275009/10.S or equal</p> <p>State number of drapes/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
Line Items #21-22 are Disposables or Accessories to purchase to go with the device being bid on line #1 ONLY. Therefore, all items must be compatible with line #1.					
21	<p>Disposable, Sterile Canister with gel (must be a min. three hundred (300) mL canister as stated in line #1: Individually wrapped sterile packaging, for one-time use, to include charcoal filters (to reduce odor). Hydrophobic construction (to prevent leaks). Includes gel (to help solidify wound exudate).</p> <p>Included in each case: *One canister (sterile fluid path) with gel, tubing, clamp and connector Packaged: 10/case</p> <p>KCI V.A.C. Freedom 300 mL Canister w/Gel, Item #320058/10.S or equal State canister size: _____</p> <p>State number of canisters/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
22	<p>Disposable, Sterile Canister with gel (must be a minimum 300 mL canister as stated in line #1: Individually wrapped sterile packaging, for one-time use, to include charcoal filters (to reduce odor). Hydrophobic construction (to prevent leaks). Includes gel (to help solidify wound exudate).</p> <p>Included in each case: *One canister (sterile fluid path) with gel, tubing, clamp and connector Packaged: 5/case</p> <p>KCI V.A.C. Freedom 300 mL Canister w/Gel, Item #320058/5.S or equal State canister size: _____</p> <p>State number of canisters/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
Line Items #23-31 are Disposables or Accessories to purchase to go with line #2 ONLY. Therefore, these items must be compatible with Line Item #2:					
23	<p>Disposable, Sterile Dressing Kit with advanced drape, Size Small (to go with line #2): Individually wrapped sterile kit, for one-time use,</p>				

Attachment B: Pricing Sheet
for Solicitation # 3000023542
VFMC Wound Vac Rental and Supplies

Line #	Item Description	Estimated Quantity	Unit of Measure	Unit Price	Extended Amount
	<p>Included in each kit:</p> <ul style="list-style-type: none"> *Two foam dressings (dimensions vary according to kit size; can be cut to size) *Two or five sheets of advanced drapes (approx. size: 8.25"X10.25", number of sheets depends on kit size; can be cut to size) *One foam quantity label *Two or four packets of no sting barrier (number of packets depends on kit size) *One pad and tubing set *One disposable ruler <p>The advanced drape to have stronger adhesive properties than the standard drape. Designed to maintain a seal even in moisture-rich environments. Dressing size: approx. 11 X 8 X 1.8 cm Contains at least two sheets of advanced drape Contains at least two packets of barrier Packaged: 5 kits/case</p> <p>KCI V.A.C. VeraFlo Dressing, Small, Item #ULTVFL05SM.S or equal</p> <p>State number of kits/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
24	<p>Disposable, Sterile Dressing Kit with advanced drape, Size Medium (to go with line #2): Individually wrapped sterile kit, for one-time use, Included in each kit:</p> <ul style="list-style-type: none"> *Two foam dressings (dimensions vary according to kit size; can be cut to size) *Two or five sheets of advanced drapes (approx. size: 8.25"X10.25", number of sheets depends on kit size; can be cut to size) *One foam quantity label *Two or four packets of no sting barrier (number of packets depends on kit size) *One pad and tubing set *One disposable ruler <p>The advanced drape to have stronger adhesive properties than the standard drape. Designed to maintain a seal even in moisture-rich environments. Dressing size: approx. 17 X 15 X 1.8 cm Contains at least two sheets of advanced drape Contains at least two packets of barrier Packaged: 5 kits/case</p> <p>KCI V.A.C. VeraFlo Dressing, Medium, Item #ULTVFL05MD.S or equal</p> <p>State number of kits/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
25	<p>Disposable, Sterile Dressing Kit with Advanced Drape, Size Large (to go with line #2): Individually wrapped sterile kit, for one-time use, Included in each kit:</p> <ul style="list-style-type: none"> *Two foam dressings (dimensions vary according to kit size; can be cut to size) *Two or five sheets of advanced drapes (approx. size: 8.25"X10.25", number of sheets depends on kit size; can be cut to size) *One foam quantity label *Two or four packets of no sting barrier (number of packets depends on kit size) *One pad and tubing set *One disposable ruler <p>The advanced drape to have stronger adhesive properties than the standard drape. Designed to maintain a seal even in moisture-rich environments. Dressing size: approx. 25 X 15 X 1.6 cm Contains at least five sheets of advanced drape Contains at least four packets of barrier Packaged: 5 kits/case</p> <p>KCI V.A.C. VeraFlo Dressing, Large, Item #ULTVFL05LG.S or equal</p> <p>State number of kits/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
26	<p>Disposable, Sterile Cleanse Dressing Kit (to go with line #2): Individually wrapped sterile kit, for one-time use. Size medium To include:</p>				

Attachment B: Pricing Sheet
for Solicitation # 3000023542
VFMC Wound Vac Rental and Supplies

Line #	Item Description	Estimated Quantity	Unit of Measure	Unit Price	Extended Amount
	<p>*One foam dressing (approx. 61 X 3.2 cm diameter) *Two sheets of advanced drapes (approx. size: 8.25"X10.25", can be cut to size) *One foam quantity label *Two packets of no sting barrier *One pad and tubing set *One disposable ruler Packaged: 5 kits/case KCI VeraFlo Cleanse Dressing, Medium, Item #ULTVCL05MD.S or equal State number of kits/case: _____ Specify brand (and number if applicable): _____</p>	1	Case		
27	<p>Disposable, Sterile Cassette (to go with line #2): Individually wrapped, for one-time use. To include instillation tubing. Connects the therapy unit to the dressing tubing and to the user-provided solution bag/bottle. Provides solution storage and delivery. Packaged: 5/case KCI V.A.C. VeraLink Cassette, Item #ULTLNK0500.S or equal State number of cassettes/case: _____ Specify brand (and number if applicable): _____</p>	1	Case		
28	<p>Disposable, Sterile Duo Tubing Set (to go with line #2): Individually wrapped, for one-time use. Separate delivery and removal pads to help enhance fluid circulation . Packaged: 5 tubing sets/case KCI V.A.C. VeraT.R.A.C. Duo Tube Set, Item #ULTDUO0500.S or equal State number of tubing kits/case: _____ Specify brand (and number if applicable): _____</p>	1	Case		
29	<p>Disposable, Sterile Canister with Gel Must be a min. 500 mL Canister as stated in line #2 Individually wrapped sterile packaging, for one-time use. To include charcoal filters (to reduce odor). Hydrophobic construction (to prevent leaks). Includes gel (to help solidify wound exudate). Includes one canister (sterile fluid path) with gel, tubing, clamp and connector . Packaged: 5/case KCI InfoV.A.C. 500mL Canister w/Gel, Item #M8275063/5.S or equal State canister size: _____ State number of canisters/case: _____ Specify brand (and number if applicable): _____</p>	1	Case		
30	<p>Disposable, Sterile Canister with Gel Must be a min. 500mL Canister as stated in line # 2 Individually wrapped sterile packaging, for one-time use. To include charcoal filters (to reduce odor). Hydrophobic construction (to prevent leaks). Includes gel (to help solidify wound exudate). Includes one canister (sterile fluid path) with gel, tubing, clamp and connector. Packaged: 10/case KCI InfoV.A.C. 500mL Canister w/Gel, Item #M8275063/10.S or equal State size of canister: _____ State number of canisters/case: _____</p>	1	Case		

Attachment B: Pricing Sheet
for Solicitation # 3000023542
VFMC Wound Vac Rental and Supplies

Line #	Item Description	Estimated Quantity	Unit of Measure	Unit Price	Extended Amount
	Specify brand (and number if applicable): _____				
31	<p>Disposable, Sterile canister with Gel Must be a minimum 1000 mL canister as stated in line #2</p> <p>Individually wrapped sterile packaging, for one-time use. To include charcoal filters (to reduce odor). Hydrophobic construction (to prevent leaks). Includes gel (to help solidify wound exudate). Includes one canister (sterile fluid path) with gel, tubing, clamp and connector . Packaged: 5/case</p> <p>KCI InfoV.A.C. 1000mL Canister, Item #M8275093/5.S or equal</p> <p>State size of canister: _____</p> <p>State number of canisters/case: _____</p> <p>Specify brand (and number if applicable): _____</p>	1	Case		
32	<p>Purchase of a Disposable, Sterile Incision Management System: Specifications: Individually wrapped sterile packaging, for one-time use.</p> <p>Battery powered, single-patient use, lightweight, portable, one-touch operation (pre-set at -125 mmHg), approx. 45 mL canister incisions management system.</p> <p>Intended to manage the environment of surgical incisions that continue to drain following sutured or stapled closure by maintaining a closed environment and removing exudate via the application of negative pressure wound therapy (NPWT) (holds incision together, protects the incision site from external infectious sources, and removes fluid and infectious materials from the surgical site).</p> <p>Can be applied for a maximum of seven days.</p> <p>Automatic time-out after 192 hours (eight days of cumulative run time).</p> <p>Can be used with both drain tubes and pain devices.</p> <p>Packaged as an all-in-one system (nothing else is needed).</p> <p>Included in each System:</p> <ul style="list-style-type: none"> - Integrated one-piece "peel & place" self-adhesive dressing with built-in pressure indicator (an adhesive dressing for use over the surgical incision and surrounding intact skin, designed to conform to patient contours: cutting/sizing not typically required). - Removable approx. 45mL canister (a sterile reservoir for collection of wound fluids). - Battery powered therapy unit (delivers negative pressure to the surgical area), including batteries pre-installed. - Carrying case (to facilitate patient mobility) with adjustable strap (for wearing over shoulder or across chest). Can also be worn on belt. - Connector (used to connect dressing to a therapy unit). - Customizable folded dressing (specially designed dressing for application to the surgical area which contains a skin interface layer that includes 0.01% ionic silver; can be used on linear incisions, non-linear incisions and intersecting incisions), approx. 10 X 3 X 1.75 cm, will cover an incision up to 8" in length. - Patch strips (used to help seal leaks by creating a continuous adhesive seal around the customizable dressing). - One drape (used to cover dressing and provide a seal for negative pressure; can be cut to size). - Disposable ruler with removable label and oval cutout (ruler is used to measure incision before cutting; label can be used to record date of dressing application or removal; oval cutout is used as a template for cutting a hole in the drape for the interface pad). - Interface pad (used to connect the customizable dressing to the therapy unit). <p>To have Indicators and Alerts such as, but not limited to: Caution light, low battery light, audio pause light, on/off and audio pause button.</p> <p>Visual alerts with solid LEDs that cannot be turned off by the user; will only stop when the alert condition has been corrected. Audible alerts with repeated beeps that can increase in volume and can be temporarily muted and will re-occur after approx. 60 minutes unless alert condition has been corrected.</p>	1	Each		

Attachment B: Pricing Sheet
for Solicitation # 3000023542
VFMC Wound Vac Rental and Supplies

Line #	Item Description	Estimated Quantity	Unit of Measure	Unit Price	Extended Amount
	<p>The following Audible Visual Alerts must be included, but are not limited to:</p> <p>Leak alert, canister full alert, low battery alert, system error alert, device lifecycle expired</p> <p>Leak Alert must calculate the leak rate and display the leak rate of the system every few seconds: must have best, good, and marginal leak rate statuses.</p> <p>Storage Conditions: Temperature range: -4 degrees F to 140 degrees F</p> <p>Operating Conditions: Temperature range: 41 degree F to 104 degree F</p> <p>Classifications: Type BF Applied Part IEC 60601-1:2005 (3rd edition) IP24</p> <p>State a contact number for ordering: _____</p> <p>Brand/Model: KCI Prevena Peel & Place System, Item #PRE1001US.S (Single Unit) or equal</p> <p>Specify brand (and number if applicable): _____</p>				